





	loteili	expres	ssion .				
	Possible yield (mg)						
Expression	mg/ml	50 ml culture	100 ml culture	250 ml culture	500 ml culture	1000 ml culture	
low	0.25	0.125	0.25	0.63	1.25	2.5	
mid	0.4	0.40	0.80	2.0	4.0	8.0	
mid-strong	0.7	0.70	1.4	3.5	7.0	14.0	
strong	1	1.0	2.0	5.0	10.0	20.0	
very strong	5	5.0	10	25	50	100	



















An	tibiotic	Stock Concentration	Storage	Working Conc (dilution)	
An	npicillin (Sodium Salt)	50 mg/ml in water (500X)	-20°C	100µg/ml (2 µl of stock/ml)	
Ch	loramphenicol	34 mg/ml in EtOH (200X)	-20°C	170 μg/ml (5 μl of stock/ml)	
Ka	namycin	25 mg/ml in water (500x)	-20°C	50 μg/ml (2 μl of stock/ml)	
Str	eptomycin	10 mg/ml in water (200X)	-20°C	50 μg/ml (5 μl of stock/ml)	
Te	tracycline HCl	5 mg/ml in EtOH (100X)	-20°C	50 µg/ml (10 µl of stock/ml)	
To maintain of membranes of drug works by (usually a hyci antibiotic in t Culturing Cel you can use fi yolu can use fo), plates are old NOT let the fi	tells that only carry you of <i>E</i> . <i>coli</i> by inhibiting th y blocking protein synti irolysis of the compour he freezer. While it is i s – Starting your cultur rom glycerol stock but . They may look good i l, either transform a ne rozen cell thaw) and so	IT plasmid, an antibiotic should ne crossilinking of the bacterial hesis at the mRNA level. It is is nd) above 60oC or if left at roo not proper to re-freeze, we fin re should almost always be do you can lose your plasmid this but are likely dead or contamir w set of cells or chip of a bit o read on an L8 Agar plate with	I be included in al membrane. Anot mportant to reme m temp for sever d little problem r ne from an isolat way. Don't trust hated with a molo f frozen glycerol s antibiotic. Cultur	I phases. Ampicillin acts to damage the her commonly used antibiotic is kanan ember that the antibiotic will "break d a days. We typically keep concentrate e-freezing nunsed antibiotic. ad colony from a freshly streaked plate the plates that have been around for to fungus or some other nasty critter. tock from the top of the tube with a p overnight in the 37°C incubator and:	e nycin. Th wm" ed . In a pine too long (If the ipet tip (d store the

Expression

- Start with clone (resistance?) from plate or frozen glycerol stock
 - Pick a colony and place in starter culture

 - Incubate overnight (o/n) OR limit OD
 OD Optical Density diffusion of light through media (not absorbanc)
 - + Some cells stop expressing well if OD is > 0.7 units











pansion of Starter Culture

A state is during the adding the starter culture larger volume LB with antibiotic (room temp) and incubate for 1-4 hours until culture density reaches 0.5-06 ODB00. Top the flask with cotton or a culture flask cap to allow culture to oxygenate without getting environmental contamination. Do NOT seal with parafilm or aluminum foil.

OPTION 1 ROOM TEMP (20°C) INDUCTION. Induce Expression – After culture has reached OD 0.5-0.6 at 600 nm, cool down culture to room temperature by placing in ice. - Induce expression by adding IPTS to a final concentration of 0.1 to 1.0 mM. IPTG is a frozen solution in the -20oC freezer. - Induce overnight (12-18 hours) at room (20oC) temp. with shaking. Refereze unused IPTG. • Note – for very difficult expression proteins you may need to express at 16°C. In that case, induce for 24-48 hours.

OPTION 2 37°C INDUCTION

- Induce Expression (see note below) After culture has reached OD 0.5-0.6 induce expression by adding IPTG to a final concentration of 0.5 mM. IPTG is a frazen solution in the -200C freezer.
 Induce for 3 to 4 hours at 37oC with shaking.
- Induce for 3 to 4 nours
 Refreeze unused IPTG.





Induction

- Allow cells to grow to moderate density before inducing e. coli to express protein
- 0.4-0.6 OD before starting culture
- Cells grow exponentially! Growth starts slowly and speeds up fast at the end. DON'T overgrow! Impacts health of cells and expression levels.
- Some plasmids are leaky may need to add glucose to decrease expression until IPTG is added



Cell Harvesting

Carefully balance tubes and centrifuge cells for pellet – can be frozen for later processing







Cell Lysis

- Low equipment/novel approach autolysis • XJ autolysis cells (BL21 strain) from Zymo Research - with second arabinose induced plasmid driving expression of bacteriophage gamma endolysin
 - Endolysin hydrolytic enzyme degrading cell wall
 Allows cell to be damage after freeze-thaw
- Lysozyme/Sonnication (without TX-100)
 - Addition of PMSF to inhibit released proteases
 - βMercaptoethanol help maintain redox state of protein
 - DNAse A and sonication reduce DNA viscosity large preps
 - Centrifuge membrane from soluble portion = lysate













His-Tag Nickle Protein Purification

Purification Instructions (per 1000 ml culture adjust to the appropriate vol

Preparation of VM-Japanes Beack/Retim: • Prepare 28 on flowshy transferring 50 ml of a 50% slurry of beack equilibrate into a clean column. Wash and equilibr the column by running 200 ml of His Elution Buffer followed by 500 ml of His Binding Buffer through the column. This SHOULD be done ahead of time! Store prepared beack with a few ml of His Binding Buffer at RT.

Purification: o Save 100 <u>ul</u> of lysate. Add clarified lysate (if frozen, check for ppt material. If there is any clumpy or ppt m lysate is cloudy, centrifuge and keep the supernatant) to the washed beads.

Barch Binding - If protein expression is low or the His-tagged protein binds with low affinity then use a batch purification methor. Combine the washed beads and lystate onto a drained and capped columor or a 50 m fals on tube for smaller volumes. Use a spatula and/or a transfer pipette to suspend the beads. Tightly cover with parafilm and incubat rocker for 30 min at room temp.

Replace the column on the stand and allow most of the beads settle then open column.
 Optional: Reaply the flow thru. This is the non-binding protein. Continue with purification

Column Binding - Flow the clarified lysate through the column. Save the eluate as flow-through in one fraction. If there is a concern with binding efficiency the flow-through can be reaching of Continue with the surfication.











Post Purification Storage & Preparation

Prevent Aggregation –

 MDH @ >1mg/ml will PPT

 Dialysis for imidazole & salt removal and buffer exchange
 Glycerol and other additives for stable storage
 DO NOT FREEZE MDH without additions!!!!





