



# Protocol:

## *Bacterial Lysis for Recombinant Protein Purification*

### Cell lysate preparation.

Resuspend Pellet Volumes: 20 - 100 ml culture - resuspend in 2-10 ml; 100-500 ml 10-25 ml; 1-2 liter 25-50 ml

- Resuspend in the appropriate binding buffer (GST or His purification) with the following additions
  - o for >500 ml cultures. Add 100 ul of 5 mg/ml DNaseA (DnaseA 0.003 mg per 1 liter pellet or 2 ul of 25000U/ml per ml). If using Pierce DNase I cat 89836, use 5 ul for a liter pellet. May need to add 5mM MgCl<sub>2</sub> and 130 uM Ca<sup>+2</sup> for optimal activity.
    - Optional 100X DNase I reaction buffer 100 mM Tris pH 7.5, 250 mM MgCl<sub>2</sub>, 10 mM CaCl<sub>2</sub>, 500 mM EDTA.
  - o Add EDTA free protease inhibitor as per manufacture's instruction. Once dissolved, the half life is quick so do not add until ready. OR 1 mM PMSF depending on material access
  - o 10 mM β (2-) mercaptoethanol \*\*\* Do not add until the day of use. (*pure β ME is 14 M*)
  - o Depending on the pellet viscosity you may need to sonicate to get the pellet fully suspended.

AutoLysis Instructions: Resuspend pellet up and down with a transfer pipette until the mixture is homogeneous

XJ Autolysis™ *E. coli* strains are a new alternative for bacterial transformation and lysis. These strains are efficiently lysed following arabinose-induced expression of the bacteriophage λ endolysin protein, coupled to a single freeze-thaw cycle.

*For optimal lysis, more care needs to be taken when selecting the lysis buffer. However, even very low concentrations of a detergent improve lysis significantly XJb is ideal for recombinant protein expression.*

- Freeze on ice or isopropanol/dry ice bath. Thaw, resuspend and freeze a second time.
- Fully resuspend pellet with disposable pipettor.
- Centrifuge at 8,000 x g for 20 min at 4°C.
- Transfer supernatant to clean 50 ml falcon tube(s). Continue purification.

### Standard, Non - Autolysis Instructions

- o Add Lysozyme (final conc of 1.0-0.5 mg/ml] - Refreeze unused lysozyme).
- o Optional - 0.1% Triton X-100 (**for prevention of aggregation of hydrophobic and membrane proteins**). **Detergents chosen for the lysis solution should be specific to the proteins**
- o Incubate on ice for 30 min while rocking / resuspending manually.
- o Sonicate ON ICE for 3 x 1-3 min bursts as high as your sample can take without cavitation.
- o Centrifuge at 8,000 x g for 20 min at 4°C. NOT in falcon tubes. Use the smaller polypropylene tubes. Must be filled at least halfway to centrifuge at this speed.
- o Transfer supernatant to clean 50 ml falcon tube(s). Continue purification.