

Single Cell Cloning by Serial Dilution



*John A. Ryan, Ph.D.
Corning Incorporated
Life Sciences
45 Nagog Park
Acton, MA 01720*

Introduction

This technique is widely used for clonal isolation of hybridomas and other cell lines that are not attachment dependent. However, it is also very useful for cloning attachment dependent cells when the cell plating efficiency is very low, unknown or unpredictable. This method is fast and easy; however, like most clonal isolation methods, there is no guarantee that the colonies arose from single cells. Recloning a second time is advised to increase the likelihood that the cells originated from a single cell.

Supplies

Nonsterile

1. Pipetting aids - Corning Cat. # 4910 (1)
2. Disposal tray or bucket for used pipettes (1)
3. Marking pen (1)
4. 200 μ L pipettor - Corning Cat. # 4963 (1)
5. 8-channel multichannel 200 μ L pipettor- Corning Cat. # 4888 (1)

Sterile

1. Cell culture medium – (Appropriate culture medium for the cells that will be cloned) (30mL)
2. Cell suspension at 5×10^4 to 1×10^5 cells/mL (200 μ L)
3. 96 well cell culture plate - Corning Cat. # 3585 (1)




Corning offers a variety of multichannel pipettors for use in 96 well plates.

4. Sterile pipettor tips
5. Reagent dispensing reservoir/tray - Corning Cat. # 4870 or 4871 (1)
6. 1, 5, and 10mL pipettes - Corning Cat. # 4485, 4487 and 4488

Procedure

1. Fill the reagent dispensing tray with 12mL of the appropriate culture medium, then using the 8 channel pipettor add 100 μ L medium to all the wells in the 96 well plate except well A1 (see diagram below) which is left empty.

Figure 1. Initial plate setup.



Adding 1000 to 2000 cells in well A1 (5×10^3 to 1×10^4 cells/mL) is a good starting cell concentration. Increase this concentration for more difficult to grow cell lines.

2. Add 200 μ L of the cell suspension to well A1. (See Figure 1.) Then using the single channel pipettor quickly transfer 100 μ L from the first well to well B1 and mix by gently pipetting. Avoid bubbles. Using the same tip, repeat these 1:2 dilutions down the entire column, discarding 100 μ L from H1 so that it ends up with the same volume as the wells above it.
3. With the 8-channel pipettor add an additional 100 μ L of medium to each well in column 1 (giving a final volume of cells and medium of 200 μ L/well). Then using the same pipettor quickly transfer 100 μ L from the wells in the first column (A1 through H1) to those in the second column (A2 through H2) and mix by gently pipetting. Avoid bubbles.
4. Using the same tips, repeat these 1:2 dilutions across the entire plate, discarding 100 μ L from each of the wells in the last column (A12 through H12) so that all the wells end up with 100 μ L of cell suspension.
5. Bring the final volume of all the wells to 200 μ L by adding 100 μ L medium to each well. Then label the plate with the date and cell type. Incubate plate undisturbed at 37°C in a humidified CO₂ incubator.

Adding filtered conditioned medium (medium in which cells have been previously grown) to the wells can increase the cloning efficiency for difficult to grow cells.

Transferring clones directly from a well in a 96 well plate into a T-25 flask is not recommended. The cells may be unable to grow due to their inability to condition the larger volume of medium in the flask. Using some conditioned medium when subculturing the cells for the first time will also help them survive and grow.

- Clones should be detectable by microscopy after 4 to 5 days and be ready to score after 7 to 10 days, depending on the growth rate of the cells. Check each well and mark all wells that contain just a single colony. These colonies can then be subcultured from the wells into larger vessels. Usually each clone is transferred into a single well in a 12 well or 24 well plate.



Figure 2. Example of a stained (1% crystal violet) plate containing a dilution plating of CHO-K1 cells. The highest cell densities occur in the wells immediately surrounding the A1 position. Wells A10, D10, E6, E11, F7, G6 and H4 appear to contain single colonies.

For additional product or technical information, please visit our web site at www.corning.com/lifesciences or call at 1-800-492-1110. International customers can call at 978-635-2200.

Acknowledgements

This protocol has evolved from a protocol I developed for cell culture training courses at the University of Connecticut, Storrs, Connecticut. I would like to thank all of my colleagues and students who, over the years, have contributed ideas and suggestions to its development.

Corning Incorporated Life Sciences

45 Nagog Park
Acton, MA 01720
t 800.492.1110
t 978.635.2200
f 978.635.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA

Australia
t 61 2-9416-0492
f 61 2-9416-0493

China
t 86 21-3222-4666
f 86 21-6288-1575

Hong Kong
t 852-2807-2723
f 852-2807-215

India
t 91 11 341 3440
f 91 11 341 1520

Japan
t 81 (0) 3-3586 1996
f 81 (0) 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6735-2913

Taiwan
t 886 2-2716-0338
f 886 2-2716-0339

EUROPE

France
t 0800 916 882
f 31 20 659 7673

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 (0) 20 655 79 28
f 31 (0) 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

LATIN AMERICA

Brazil
t (55-11) 3089-7420
f (55-11) 3167-0700

Mexico
t (52-81) 8313-8586
f (52-81) 8313-8589