



pcDNA[™] 3.1/Hygro (+)

pcDNA[™] 3.1/Hygro (-)

**For high-level stable and transient expression in
mammalian hosts**

Catalog nos. V870-20, V875-20

Version G

10 November 2010

28-0149

User Manual

Table of Contents

Kit Contents and Storage	v
Accessory Products.....	vi
Introduction	1
Overview	1
Methods	2
Cloning into pcDNA™3.1/Hygro (+/-).....	2
Transfection	6
Creation of Stable Cell Lines	7
Appendix	10
Map of pcDNA™3.1/Hygro (+) and pcDNA™3.1/Hygro (-) Vectors	10
Features of pcDNA™3.1/Hygro (+) and pcDNA™3.1/Hygro (-) Vectors.....	11
Map of pcDNA™3.1/Hygro/ <i>lacZ</i>	12
Technical Support	13
Purchaser Notification	14
References	16

Kit Contents and Storage

Shipping and Storage

pcDNA[™]3.1/Hygro (+/-) vectors are shipped on wet ice. Upon receipt, store vectors at -20°C.

Kit Contents

Each catalog number contains the following vectors. All vectors are supplied in aliquot detailed below. **Store the vectors at -20°C.**

Catalog nos.	Vector	Quantity	Composition (supplied as)
V870-20	pcDNA [™] 3.1/Hygro (+)	20 µg	40 µl of 0.5 µg/µl pcDNA [™] 3.1/Hygro (+) vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0
	pDNA [™] 3.1/Hygro/ <i>lacZ</i> control	20 µg	40 µl of 0.5 µg/µl pcDNA [™] 3.1/Hygro/ <i>lacZ</i> control vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0
V875-20	pcDNA [™] 3.1/Hygro (-)	20 µg,	40 µl of 0.5 µg/µl pcDNA [™] 3.1/Hygro (-) vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0
	pDNA [™] 3.1/Hygro/ <i>lacZ</i> control	20 µg	40 µl of 0.5 µg/µl pcDNA [™] 3.1/Hygro/ <i>lacZ</i> control vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0

Accessory Products

Introduction

The following additional products may be used with the pcDNA™3.1/Hygro (+/-) vectors. For more information, visit www.invitrogen.com or contact **Technical Support** (see page 13).

Item	Quantity	Catalog no.
One Shot® TOP10F' (chemically competent cells)	21 × 50 µl	C3030-03
One Shot® TOP10 Chemically Competent <i>E. coli</i>	10 reactions	C4040-10
One Shot® TOP10 Electrocomp™ <i>E. coli</i>	10 reactions 20 reactions	C4040-50 C4040-52
One Shot® MAX Efficiency® DH10B™ (chemically competent cells)	1 ml	18297-010
One Shot® MAX Efficiency® DH10α™ (T1 _R competent cells)	20 × 50 µl	12297-016
T7 Promoter Primer	2 µg	N560-02
BGH Reverse Primer	2 µg	N575-02
PureLink™ HiPure Plasmid Miniprep Kit	100 preps	K2100-03
PureLink™ HiPure Plasmid Midiprep Kit	25 preps	K2100-04
Hygromycin B	20 ml	10687-010
<i>Bgl</i> II	400 units	15213-010
Ssp	400 units	15458-011
β-Gal Assay Kit	1 kit	K1455-01
β-Gal Staining Kit	1 kit	K1465-01

Introduction

Overview

Introduction

pcDNA™3.1/Hygro (+) and pcDNA™3.1/Hygro (-) are 5.6 kb vectors derived from pcDNA™3.1 and are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. These vectors contain the following elements:

- Human cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells
- Multiple cloning sites in the forward (+) and reverse (-) orientations to facilitate cloning
- Hygromycin resistance gene for selection of stable cell lines
- Episomal replication in cells lines that are latently infected with SV40 or that express the SV40 large T antigen (e.g., COS-1, COS-7)

A control plasmid, pcDNA™3.1/Hygro/*lacZ*, is included for use as a positive control for transfection and expression in the cell line of choice.

CMV Promoter

pcDNA™3.1/ Hygro (+) and pcDNA™3.1/ Hygro (-) vectors contain the human CMV immediate early promoter to allow high-level, constitutive expression of the gene of interest in mammalian cells (Andersson *et al.*, 1989; Boshart *et al.*, 1985; Nelson *et al.*, 1987). Although highly active in most mammalian cell lines, activity of the viral promoter can be down-regulated in some cell lines due to methylation (Curradi *et al.*, 2002), histone deacetylation (Rietveld *et al.*, 2002), or both.

Experimental Outline

Use the following outline to clone and express your gene of interest in pcDNA™3.1/Hygro (+/-).

1. Consult the multiple cloning sites described on pages 3–4 to design a strategy to clone your gene into pcDNA™3.1/Hygro.
 2. Ligate your insert into the appropriate vector and transform into *E. coli*. Select transformants on LB plates containing 50–100 µg/ml ampicillin.
 3. Analyze your transformants for the presence of insert by restriction digestion.
 4. Select a transformant with the correct restriction pattern and sequence to confirm that your gene is cloned in the proper orientation.
 5. Transfect your construct into the mammalian cell line of interest using your own method of choice. Generate a stable cell line, if desired.
 6. Test for expression of your recombinant gene by western blot analysis or functional assay.
-

Methods

Cloning into pcDNA™ 3.1/Hygro (+/–)

Introduction

To recombine your gene of interest into pcDNA™3.1/Hygro (+/–), you will need to ligate your gene of interest into pcDNA™3.1/Hygro (+) or pcDNA™3.1/Hygro (–). Diagrams of the multiple cloning sites for each vector are provided on pages 3–4. General considerations for cloning and transformation are listed below.

General Molecular Biology Techniques

For help with DNA ligations, *E. coli* transformations, restriction enzyme analysis, purification of single-stranded DNA, DNA sequencing, and DNA biochemistry, refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

E. coli Strain

Many *E. coli* strains are suitable for the propagation of pcDNA™3.1/Hygro (+/–) including TOP10F′, DH5α, and TOP10 (see page vi for ordering information). We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*) and endonuclease A-deficient (*endA*). For your convenience, TOP10F′ is available as chemically competent or electrocompetent cells from Invitrogen.

Transformation Method

You may use any method of your choice for transformation. Chemical transformation is the most convenient for most researchers. Electroporation is the most efficient and the method of choice for large plasmids.

Maintenance of pcDNA™ 3.1/Hygro (+/–)

To propagate and maintain pcDNA™3.1/Hygro (+/–), we recommend using 10 ng of the vector to transform a *recA*, *endA* *E. coli* strain such as TOP10, TOP10F′, DH5α™, or equivalent (see page vi for ordering information) using your method of choice. Select transformants on LB plates containing 50–100 µg/ml ampicillin.

For long-term storage, be sure to prepare a glycerol stock (page 5) of your plasmid-containing *E. coli* strain.

Points to Consider Before Recombining into pcDNA™ 3.1/Hygro (+/–)

pcDNA™3.1/Hygro (+) and pcDNA™3.1/Hygro (–) are nonfusion vectors. Your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak, 1990; Kozak, 1991). An example of a Kozak consensus sequence is provided below. Other sequences are possible, but the G or A at position –3 and the G at position +4 (shown in bold) illustrates the most commonly occurring sequence with strong consensus. Replacing one of the two bases at these positions provides moderate consensus, while having neither results in weak consensus. The ATG initiation codon is shown underlined.

(G/A)NNATGG

Your insert must also contain a stop codon for proper termination of your gene. Note that the *Xba* I site contains an internal stop codon (TCTAGA).

Continued on next page

Cloning into pcDNA™ 3.1/Hygro, Continued

Multiple Cloning Site of pcDNA™ 3.1/Hygro (+)

Below is the multiple cloning site for pcDNA™ 3.1/Hygro (+). Restriction sites are labeled to indicate the cleavage site. The *Xba* I site contains an internal stop codon (TCTAGA). The multiple cloning site has been confirmed by sequencing and functional testing. The complete sequence of pcDNA™ 3.1/Hygro (+) is available for downloading from our web site at www.invitrogen.com or by from **Technical Support** (see page 13).

```

          enhancer region (3' end)
689  CATTGACGTC AATGGGAGTT TGTTTTGGCA CAAAATCAA CGGGACTTTC CAAAATGTGC

          CAAT
749  TAACAAC TCC GCCCCATTGA CGCAAATGGG CGGTAGGCGT GTACGGTGGG AGGTCTATAT
          3' end of hCMV
          putative transcriptional start
809  AAGCAGAGCT CTCTGGCTAA CTAGAGAACC CACTGCTTAC TGGCTTATCG AAATTAATAC

T7 promoter/primer binding site
869  GACTCACTAT AGGGAGACCC AAGCTGGCTA GCGTTTAAAC TTAAGCTTGG TACCGAGCTC
          BamH I
          BstX I*
          EcoR V
          BstX I*
          Not I
          Xho I
929  GGATCCACTA GTCCAGTGTG GTGGAATTCT GCAGATATCC AGCACAGTGG CGGCCGCTCG

          Xba I
          Apa I
          Pme I
          BGH reverse priming site
989  AGTCTAGAGG GCCCGTTTAA ACCCGCTGAT CAGCCTCGAC TGTGCCTTCT AGTTGCCAGC

1049  CATCTGTTGT TTGCCCTCC CCCGTGCCTT CTTGACCCT GGAAGGTGCC ACTCCCCTG

          BGH poly (A) site
1109  TCCTTTCCTA ATAAAATGAG GAAATTGCAT
  
```

Continued on next page

Cloning into pcDNA™ 3.1/Hygro, Continued

Multiple Cloning Site of pcDNA™ 3.1/Hygro (-)

Below is the multiple cloning site for pcDNA™ 3.1/Hygro (-). Restriction sites are labeled to indicate the cleavage site. The *Xba* I site contains an internal stop codon (TCTAGA). The multiple cloning site has been confirmed by sequencing and functional testing. The complete sequence of pcDNA™ 3.1/Hygro (-) is available for downloading from our web site at www.invitrogen.com or from **Technical Support** (see page 13).

```

          enhancer region (3' end)
689  CATTGACGTC AATGGGAGTT TGTTTTGGCA CAAAATCAA CGGGACTTTC CAAAATGTCTG

          CAAT
749  TAACAAC TCC GCCCCATTGA CGCAAATGGG CGGTAGGCGT GTACGGTGGG AGGCTCTATAT
          3' end of hCMV
          putative transcriptional start
809  AAGCAGAGCT CTCTGGCTAA CTAGAGAACC CACTGCTTAC TGGCTTATCG AAATTAATAC

T7 promoter/primer binding site      Nhe I      Pme I      Apa I Xba I      Xho I      Not I
869  GACTCACTAT AGGGAGACCC AAGCTGGCTA GCGTTTAAAC GGGCCCTCTA GACTCGAGCG

          BstX I*  EcoR V
929  GCCGCCACTG TGCTGGATAT CTGCAGAATT CCACCACACT GGACTAGTGG ATCCGAGCTC

          Asp718 I  Kpn I  Hind III  Afl II      Pme I
989  GGTACCAAGC TTAAGTTTAA ACCGCTGATC AGCCTCGACT GTGCCTTCTA GTTGCCAGCC

          BGH reverse priming site
1049 ATCTGTTGTT TGCCCCTCCC CCGTGCCTTC CTTGACCCTG GAAGGTGCCA CTCCCCTGT

          BGH poly (A) site
1109 CCTTTCCTAA TAAAATGAGG AAATTGCATC
  
```

Cloning into pcDNA™ 3.1/Hygro, Continued

E. coli **Transformation**

Transform your ligation mixtures into a competent *recA*, *endA* *E. coli* strain (e.g. TOP10F', DH5 α , TOP10, page vi) and select on LB plates containing 50–100 $\mu\text{g}/\text{ml}$ ampicillin. Select 10–20 clones and analyze for the presence and orientation of your insert.



We recommend that you sequence your construct with the T7 Promoter and BGH Reverse primers (page vi) to confirm that your gene is in the correct orientation for expression, and contains an ATG initiation codon and a stop codon. Refer to the diagrams on pages 3–4 for the sequences and location of the priming sites.

Primer	Sequence
BGH Reverse	5'-TAGAAGGCACAGTCGAGG-3'
T7 Promoter	5'-TAATACGACTCACTATAGGG-3'

For your convenience, Invitrogen offers a custom primer synthesis service. Visit www.invitrogen.com for more details.

Preparing a Glycerol Stock for Long-Term Storage

Once you have identified the correct clone, purify the colony and make a glycerol stock for long-term storage. You should keep a DNA stock of your plasmid at -20°C .

1. Streak the original colony out on an LB plate containing 50 $\mu\text{g}/\text{ml}$ ampicillin. Incubate the plate at 37°C overnight.
 2. Isolate a single colony and inoculate into 1–2 ml of LB containing 50 $\mu\text{g}/\text{ml}$ ampicillin.
 3. Grow the culture to mid-log phase ($\text{OD}_{600} = 0.5\text{--}0.7$).
 4. Mix 0.85 ml of culture with 0.15 ml of sterile glycerol and transfer to a cryovial.
 5. Store at -80°C .
-

Continued on next page

Transfection

Introduction

Once you have verified that your gene is cloned in the correct orientation and contains an initiation ATG and a stop codon, you are ready to transfect your cell line of choice. We recommend that you include the positive control vector and a mock transfection (negative control) to evaluate your results.

Plasmid Preparation

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells, and salt will interfere with lipids, decreasing transfection efficiency. We recommend isolating plasmid DNA using the PureLink™ HiPure Miniprep Kit or the PureLink™ HiPure Midiprep Kit (see page vi for ordering information), or CsCl gradient centrifugation.

Methods of Transfection

For established cell lines (e.g., HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. We recommend that you follow the protocol for your cell line, exactly. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

Methods for transfection include calcium phosphate (Chen & Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1989; Felgner & Ringold, 1989) and electroporation (Chu *et al.*, 1987; Shigekawa & Dower, 1988). For high efficiency transfection in a broad range of mammalian cells, we recommend using Lipofectamine™ 2000 Reagent available from Invitrogen. For more information on Lipofectamine™ 2000 and other transfection reagents available, visit our web site at www.invitrogen.com or contact **Technical Support** (page 13).

Positive Control

pcDNA™3.1/Hygro/*lacZ* is provided as a positive control vector for mammalian transfection and expression (see page 12) and may be used to optimize transfection conditions for your cell line. The gene encoding β -galactosidase (*lacZ*) is expressed in mammalian cells under the control of the CMV promoter. A successful transfection will result in β -galactosidase expression that can be easily assayed (see below).

Assay for β -Galactosidase Activity

You may assay for β -galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. Invitrogen offers the β -Gal Assay Kit and the β -Gal Staining Kit (see page vi for ordering information) for fast and easy detection of β -galactosidase expression.

Creation of Stable Cell Lines

Introduction

The pcDNA™3.1/Hygro (+) and pcDNA3.1/Hygro (-) vectors contain the hygromycin resistance gene for selection of stable cell lines using hygromycin B. We recommend that you test the sensitivity of your mammalian host cell to hygromycin B, as natural resistance varies among cell lines. General information and guidelines are provided in this section for your convenience.

Hygromycin-B Activity

Hygromycin-B (527.5 MW) is an aminocyclitol that inhibits protein synthesis by disrupting translocation and promoting mistranslation. Hygromycin-B-phosphotransferase detoxifies hygromycin-B by phosphorylation.



- Hygromycin is light sensitive. Store the liquid stock solution at 4°C protected from exposure to light.
 - Hygromycin is toxic. Do not ingest solutions containing the drug.
 - Wear gloves, a laboratory coat, and safety glasses or goggles when handling hygromycin and hygromycin-containing solutions.
-

Determining Antibiotic Sensitivity

To successfully generate a stable cell line expressing your gene of interest from pcDNA™3.1/Hygro, you need to determine the minimum concentration of hygromycin B required to kill your untransfected host cell line. Typically, concentrations ranging from 10 to 400 g/ml hygromycin are sufficient to kill most untransfected mammalian cell lines. We recommend that you test a range of concentrations (see protocol below) to ensure that you determine the minimum concentration necessary for your host cell line.

1. Plate or split a confluent plate so the cells will be approximately 25% confluent. Prepare a set of 7 plates. Allow cells to adhere overnight.
2. The next day, substitute culture medium with medium containing varying concentrations of hygromycin (0, 10, 25, 50, 100, 200, 400 g/ml hygromycin).
3. Replenish the selective media every 3–4 days, and observe the percentage of surviving cells.
4. Count the number of viable cells at regular intervals to determine the appropriate concentration of hygromycin that prevents growth within 2–3 weeks after addition of hygromycin.

Note: Cells will divide once or twice in the presence of lethal doses of hygromycin, so the effects of the drug may take several days to become apparent. Complete inhibition of cell growth can take 2–3 weeks of growth in selective medium.

Continued on next page

Creation of Stable Cell Lines, Continued

Possible Sites for Linearization of pcDNA™3.1/Hygro (+)

Prior to transfection, we recommend that you linearize the pcDNA™3.1/Hygro (+) vector. Linearizing the pcDNA™3.1/Hygro (+) will decrease the likelihood of the vector integrating into the genome in a way that disrupts the gene of interest or other elements required for expression in mammalian cells. The table below lists unique restriction sites that may be used to linearize your construct prior to transfection. **Other unique restriction sites are possible.** Be sure that your insert does not contain the restriction enzyme site you wish to use to linearize your vector.

Enzyme	Restriction Site (bp)	Location	Supplier
<i>Bgl</i> II	12	Upstream of CMV promoter	Invitrogen (page vi)
<i>Mfe</i> I	161	Upstream of CMV promoter	New England Biolabs
<i>Sap</i> I	3668	Backbone	New England Biolabs
<i>Eam</i> 1105 I	4674	Ampicillin gene	AGS*, Fermentas, Takara
<i>Ssp</i> I	5478	Backbone	Invitrogen (page vi)

*Angewandte Gentechnologie Systeme

Possible Sites for Linearization of pcDNA™3.1/Hygro (-)

The table below lists unique restriction sites that may be used to linearize your pcDNA™3.1/Hygro (-) construct prior to transfection. **Other unique restriction sites are possible.** Be sure that your insert does not contain the restriction enzyme site you wish to use to linearize your vector.

Enzyme	Restriction Site (bp)	Location	Supplier
<i>Bgl</i> II	12	Upstream of CMV promoter	Invitrogen (page vi)
<i>Mfe</i> I	161	Upstream of CMV promoter	New England Biolabs
<i>Sap</i> I	3667	Backbone	New England Biolabs
<i>Eam</i> 1105 I	4673	Ampicillin gene	AGS*, Fermentas, Takara
<i>Ssp</i> I	5477	Backbone	Invitrogen (page vi)

*Angewandte Gentechnologie Systeme

Continued on next page

Creation of Stable Cell Lines, Continued

Selection of Stable Integrants

Once you have determined the appropriate hygromycin concentration to use for selection in your host cell line, you can generate a stable cell line expressing your gene of interest.


1. Transfect your mammalian host cell line with your pcDNA™3.1/Hygro construct using the desired protocol. Remember to include a plate of untransfected cells as a negative control and the pcDNA™3.1/Hygro/*lacZ* plasmid as a positive control.
 2. 24 hours after transfection, wash the cells and add fresh medium to the cells.
 3. 48 hours after transfection, split the cells into fresh medium containing hygromycin at the pre-determined concentration required for your cell line. Split the cells such that they are no more than 25% confluent.
 4. Feed the cells with selective medium every 3–4 days until hygromycin-resistant foci can be identified.
 5. Pick and expand colonies in 96- or 48-well plates.
-

Appendix

Map of pcDNA™ 3.1/Hygro (+) and pcDNA™ 3.1/Hygro (-) Vectors

Map of pcDNA™ 3.1/Hygro (+/-)

The figure below summarizes the features of the pcDNA™ 3.1/Hygro (+) and pcDNA™ 3.1/Hygro (-) vectors. The complete sequences for pcDNA™ 3.1/Hygro (+) and pcDNA™ 3.1/Hygro (-) are available for downloading from our web site at www.invitrogen.com or from **Technical Support** (see page 13).



Comments for pcDNA3.1/Hygro (+): 5597 nucleotides

- CMV promoter: bases 209-863
- T7 promoter/priming site: bases 863-882
- Multiple cloning site: bases 895-1010
- BGH reverse priming site: bases 1022-1039
- BGH polyadenylation signal: bases 1021-1235
- f1 origin: bases 1298-1711
- SV40 promoter and origin: bases 1776-2100
- Hygromycin resistance gene: bases 2118-3141
- SV40 early polyadenylation signal: bases 3154-3526
- pUC origin: bases 3786-4456 (complementary strand)
- Ampicillin resistance gene: bases 4601-5461 (complementary strand)

Features of pcDNA™3.1/Hygro (+) and pcDNA™3.1/Hygro (-) Vectors

Features of pcDNA™3.1/Hygro (+/-)

pcDNA™3.1/Hygro (+) (5597 bp) and pcDNA™3.1/Hygro (-) (5596 bp) contain the following elements. All features have been functionally tested.


Feature	Benefit
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Permits efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert
Multiple cloning site in forward or reverse orientation	Allows insertion of your gene and facilitates cloning
BGH reverse priming site	Allows sequencing through the insert
Bovine growth hormone (BGH) polyadenylation signal	Permits efficient transcription termination and polyadenylation of mRNA (Goodwin & Rottman, 1992)
f1 origin	Allows rescue of single-stranded DNA
SV40 promoter and origin	Allows efficient, high-level expression of the hygromycin resistance gene and episomal replication in cells expressing SV40 large T antigen
Hygromycin resistance gene (Hygromycin-B-phosphotransferase)	Permits selection of stable transfectants in mammalian cells (Gritz & Davies, 1983; Palmer <i>et al.</i> , 1987)
SV40 early polyadenylation signal	Permits efficient transcription termination and polyadenylation of mRNA
pUC origin	Allows high-copy number replication and growth in <i>E. coli</i>
Ampicillin resistance gene (β -lactamase)	Permits selection of vector in <i>E. coli</i>

Map of pcDNA™ 3.1/Hygro/lacZ

Map of pcDNA™ 3.1/Hygro/lacZ

pcDNA™ 3.1/Hygro/lacZ is an 8648 bp control vector containing the gene for β-galactosidase. It was constructed by cloning a 3.2 kb *Hind* III-*Xho* I fragment containing the *lacZ* gene into pcDNA™ 3.1/Hygro (+).

The figure below summarizes the features of the pcDNA™ 3.1/Hygro/lacZ vector. The complete nucleotide sequence for pcDNA™ 3.1/Hygro/lacZ is available for downloading from our web site at www.invitrogen.com or by from **Technical Support** (see page 13).



Comments for pcDNA3.1/Hygro(+)/lacZ 8648 nucleotides

- CMV promoter: bases 209-863
- T7 promoter/priming site: bases 863-882
- LacZ ORF: bases 972-4044
- BGH reverse priming site: bases 4073-4090
- BGH polyadenylation signal: bases 4072-4286
- f1 origin: bases 4349-4762
- SV40 promoter and origin: bases 4827-5151
- Hygromycin resistance gene: bases 5169-6192
- SV40 early polyadenylation signal: bases 6205-6577
- pUC origin: bases 6837-7507 (complementary strand)
- Ampicillin resistance gene: bases 7652-8512 (complementary strand)

Technical Support

Web Resources



Visit the Invitrogen web site at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, MSDSs, FAQs, formulations, citations, handbooks, etc.
 - Complete technical support contact information
 - Access to the Invitrogen Online Catalog
 - Additional product information and special offers
-

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our web site (www.invitrogen.com).

Corporate Headquarters:

Invitrogen Corporation
5791 Van Allen Way
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail: tech_support@invitrogen.com

Japanese Headquarters:

Invitrogen Japan
LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Tel: 81 3 5730 6509
Fax: 81 3 5730 6519
E-mail: jpinfo@invitrogen.com

European Headquarters:

Invitrogen Ltd
Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail: eurotech@invitrogen.com

MSDS

MSDSs (Material Safety Data Sheets) are available on our web site at www.invitrogen.com/msds.

Certificate of Analysis

Product qualification is described in the Certificate of Analysis (CofA), available on our website by product lot number at www.invitrogen.com/cofa.

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

Invitrogen warrants that all of its products will perform according to specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore, Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

Introduction

Use of the pcDNA™3.1/Hygro (+/-) vectors is covered under a number of different licenses including those detailed below.

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008; Phone (760) 603-7200 or e-mail: outlicensing@lifetech.com.

References

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989) Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* *264*, 8222-8229
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994) *Current Protocols in Molecular Biology*, Greene Publishing Associates and Wiley-Interscience, New York
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985) A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* *41*, 521-530
- Chen, C., and Okayama, H. (1987) High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* *7*, 2745-2752
- Chu, G., Hayakawa, H., and Berg, P. (1987) Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nucleic Acids Res.* *15*, 1311-1326
- Curradi, M., Izzo, A., Badaracco, G., and Landsberger, N. (2002) Molecular Mechanisms of Gene Silencing Mediated by DNA Methylation. *Mol. Cell. Biol.* *22*, 3157-3173
- Felgner, P. L., Holm, M., and Chan, H. (1989) Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121
- Felgner, P. L. a., and Ringold, G. M. (1989) Cationic Liposome-Mediated Transfection. *Nature* *337*, 387-388
- Goodwin, E. C., and Rottman, F. M. (1992) The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. *J. Biol. Chem.* *267*, 16330-16334
- Gritz, L., and Davies, J. (1983) Plasmid-Encoded Hygromycin-B Resistance: The Sequence of Hygromycin-B-Phosphotransferase Gene and its Expression in *E. coli* and *S. Cerevisiae*. *Gene* *25*, 179-188
- Kozak, M. (1987) An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nucleic Acids Res.* *15*, 8125-8148
- Kozak, M. (1990) Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305
- Kozak, M. (1991) An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. *J. Cell Biology* *115*, 887-903
- Miller, J. H. (1972) *Experiments in Molecular Genetics*, Cold Spring Harbor Laboratory, Cold Spring Harbor, New York
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987) Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Molec. Cell. Biol.* *7*, 4125-4129
- Palmer, T. D., Hock, R. A., Osborne, W. R. A., and Miller, A. D. (1987) Efficient Retrovirus-Mediated Transfer and Expression of a Human Adenosine Deaminase Gene in Diploid Skin Fibroblasts from an Adenosine-Deficient Human. *Proc. Natl. Acad. Sci. U.S.A.* *84*, 1055-1059
- Rietveld, L. E., Caldenhoven, E., and Stunnenberg, H. G. (2002) In vivo Repression of an Erythroid-Specific Gene by Distinct Corepressor Complexes. *EMBO J.* *21*, 1389-1397
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Second Ed., Cold Spring Harbor Laboratory Press, Plainview, New York
- Shigekawa, K., and Dower, W. J. (1988) Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* *6*, 742-751
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977) Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* *11*, 223-232

©1999–2008, 2010 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.



Corporate Headquarters

Invitrogen Corporation

5791 Van Allen Way

Carlsbad, CA 92008

T: 1 760 603 7200

F: 1 760 602 6500

E: tech_support@invitrogen.com

For country-specific contact information, visit our web site at www.invitrogen.com